DICOM CD/DVD startup guide (ENGLISH)

This short guide explains how to launch the viewer application which allows displaying the DICOM medical images present on this CD/DVD. Depending on the configuration of the application which created this CD/DVD, the viewer application may also be absent from the present media.

Instructions for each supported operating system follow.

Windows operating system

If your PC runs a Windows operating system, normally the *LocalEye Viewer* application will start automatically once the CD/DVD is inserted in the drive. If this doesn't happen, please browse the CD/DVD contents and double-click on the "*LocalEyeLauncher.exe*" file located in the "*Viewer-Windows*" folder of the CD/DVD: the viewer application will start.

Please refer to the Online Help of the viewer application for further information about its usage.

Mac OS X operating system

If you are using a Mac, please insert the CD/DVD in the appropriate drive, wait a few seconds, then double click on the CD/DVD icon which will appear on your desktop, in order to open the CD/DVD and browse its contents. Double-click on

the "*ViewCDMacOSX.command*" file located in the root of the CD/DVD: the viewer application will start.

Please refer to the Online Help of the viewer application for further information about its usage.

Linux operating system

If you are using a Linux PC with a graphic windowing environment, please insert the CD/DVD in the appropriate drive, wait a few seconds, then double click on the CD/DVD icon which will appear on your desktop (auto-mounting shall be enabled), in order to open the CD/DVD and browse its contents. Double-click on the "*ViewCDLinux.sh*" file located in the root of the CD/DVD: the viewer application will start.

Please refer to the Online Help of the viewer application for further information about its usage.

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